

AD-A058 800

IMASS INC ACCORD MA
DEVELOPMENT OF ELECTRONICS, MECHANICAL FIXTURES AND COMPUTER SO--ETC(U)
AUG 78 E J TOLLE
1-78

F/G 11/4

N00014-77-C-0421

NL

UNCLASSIFIED

1 OF
AD
A058800



AD A058800

DDC FILE COPY

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER N00014-77C-0421 #1-78 ✓	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) Development of Electronics, Mechanical Fixtures and Computer Software for Automation of Dynamic Mechanical Measurements on Polymer Solids and Computer Reduction of Data		5. TYPE OF REPORT & PERIOD COVERED Technical Report - Final
7. AUTHOR(s) E. J. Tolle, President Imass, Inc.		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS Imass, Inc. ✓ P.O. Box 134 Accord (Hingham) MA 02018		8. CONTRACT OR GRANT NUMBER(s) N00014-77-C-0421 <i>ew</i>
11. CONTROLLING OFFICE NAME AND ADDRESS Office of Naval Research Arlington, Virginia		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE 28 August 1978
		13. NUMBER OF PAGES 50
		15. SECURITY CLASS. (of this report) Unclassified
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Dynamic Mechanical Modulus Tan delta ?		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) The objective of this work is to automate data taking and data reduction on a Rheovibron Dynamic Viscoelastometer. A computer program with the necessary hardware has been developed to control tension on the sample, to monitor power to heaters, to control a temperature scan, and to acquire data on polymeric solids from which temperature dependence of a dynamic modulus and tangent of the loss angle ($\tan \delta$) can be calculated. <i>Ldelta</i>		

DD FORM 1 JAN 73 1473

EDITION OF 1 NOV 65 IS OBSOLETE
S/N 0102-LF-014-6601

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

Report **N00014-77C-0421** #1-78

6 Development of Electronics, Mechanical Fixtures
and Computer Software for Automation of Dynamic
Mechanical Measurements on Polymer Solids and
Computer Reduction of Data

10 E. J. Tolle President
Imass, Inc.
P.O. Box 134
Accord (Hingham) MA 02018

11 28 Aug 1978

12 73 p.

9 Final Technical Report, Final

Approved for public release; distribution unlimited.

Office of Naval Research
Arlington, Virginia

Contract #N00014-77-C-0421

slc

393 489

mt

TABLE OF CONTENTS

	<u>Page</u>
Introduction	1
Dynamic Mechanical Measurements with the Rheovibron	2
Objectives of This Work	3
Description of the Automation System	3
Measurement Features:	
Temperature	5
Modulus and Phase Measurement	6
Sample Tension Control	7
Typical Test Procedure	8
Conclusions and Recommendations	9
Cooldown Routine	10
Main Control Routine	11
Data Reduction Routine	12
<u>Parts Lists & Diagrams</u>	
Main Frame	16
Cabinet	17
Motor Drive Assembly	18
Cables	22
Interconnecting Harness	23
Amp.-Phase-Temp. card	24
Tension card	27
Motor card	30
Range card	33
Rheovibron Relay Board	36
Rheovibron Interface	40
Heating Control	42
Power Supply Panel	45
Mounting for Temperature Chamber	47
References	50

ACCESSION for	
NTIS	Write Section <input checked="" type="checkbox"/>
DDC	B II Section <input type="checkbox"/>
UNANNOUNCED	<input type="checkbox"/>
JUSTIFICATION	
BY	
DISTRIBUTION/AVAILABILITY CODES	
In	SPECIAL
A	

INTRODUCTION

Dynamic mechanical measurements on polymers and other viscoelastic materials are an attractive means of measuring fundamental relaxations which reflect molecular motions in the sample. Such measurements are useful in studying effects of many parameters such as molecular weight, crystallinity, crosslinking, additives and plasticizers, copolymerization and blending on mechanical performance.


One of the most popular commercial instruments for making dynamic measurements is the Rheovibron which was designed and developed by Professor Motowo Takayanagi of Kyushu University, Fukuoka, Japan. His publications involving the instrument have given him recognition, worldwide, as a leading polymer scientist. The instrument was built and offered for sale commercially by Toyo Baldwin Co. Ltd. (formerly Toyo Measuring Instruments Co. Ltd.) of Tokyo, Japan.

The most serious obstacles to the use of the Rheovibron are its labor intensity (a single scan on one sample can require as much as 5 hours of undivided operator attention) and the high level of operator skill required to obtain reproduceable measurements.

This development was undertaken to automate the instrument both with regard to data acquisition and to data reduction. By employing a dedicated computer, we planned to improve both the accuracy and precision of measurement by eliminating human errors which are most prevalent in applying tension to the sample, and by eliminating operator decisions and incorporating such decisions in the computer program.

DYNAMIC MECHANICAL MEASUREMENTS WITH THE RHEOVIBRON

The instrument is designed to impose a sinusoidal strain, in tension, on the sample and to measure the sinusoidal stress response of the material with respect to amplitude and phase. From these data one can determine a complex modulus (E^*) at known frequency and temperature and a damping factor which is defined as the ratio of the quadrature and in-phase components of the complex modulus and which also represents the tangent of the phase angle which exists between the sinusoidal deformation and stress response of the sample.



$$E^* = \frac{\text{Dynamic Stress Amp.}}{\text{Dynamic Strain Amp.}}$$

$$\tan \delta = \frac{E''}{E'}$$

Figure 1.

The amplitude of the sinusoidal deformation is made small (1.0×10^{-3} cm, approximately .03 - 0.05% strain) to insure linear viscoelastic response. The frequency is maintained constant for each data point.

Data are usually taken over a range of temperature and sometimes over a limited range of frequencies. A typical range of temperature is -150°C to as far above the glass transition temperature as the operator's skill can achieve. The limit is usually determined by sag and creep of the sample under gravity and tension making further data taking impossible. The normal frequencies of the instrument are 3.5, 11, 35, 110 Hz.

OBJECTIVES OF THIS WORK

GENERAL

1. Reduce labor intensity of measurement by automating control functions and data taking.
2. Improve the precision of measurement especially with respect to $\tan \delta$.
3. Extend the range of $\tan \delta$ measurements.
4. Provide a computer program for data reduction.

DESCRIPTION OF THE AUTOMATION SYSTEM

The system retains and uses the standard electronics associated with manual operation of the Rheovibron. Transducer excitation, electromagnetic driver excitation, and reference signal are all furnished to the system by the standard instrument. The automatic control functions are sample tension and temperature. Data taking and data reduction are also automatically handled by the system.

The automation system is designed around a computer program which is implemented with a Hewlett-Packard 9825A Calculator, interfaced for control and data acquisition to a Hewlett-Packard 6940B Multiprogrammer. Commands in the program instruct the instrument through the multiprogrammer to initiate controls over the experiment.

Other program commands instruct the multiprogrammer to acquire data and return it to the calculator for manipulation and/or storage. A separate program is provided for data reduction and a line printer and/or plotter can be provided as optional equipment.

The standard phase angle measuring system has been replaced with a lock-in amplifier (Princeton Applied Research Model 2504) to simplify automation of the measurement, to improve resolution of small angles, and to extend the range of $\tan \delta$ measurements.

The simplified system is shown, schematically, in Figure 2. A Rheovibron interface, on command from the calculator, uses a signal from either the load or displacement transducer to the lock-in amplifier where it is measured. The measured result is passed through the Amp, Phase, Temp. and Relay cards to the Voltage Monitor where it is digitized and stored in the calculator memory. The two signals are alternately measured throughout the test run. The interface also passes the load transducer signal through the tension amplifier and relay card to the voltage monitor for digitizing and to the calculator

where it is compared for tension limits. If the tension requires correction, a signal from the calculator through the relay card commands the motor to increase or decrease sample tension the proper number of steps. Temperature at the sample is detected by a platinum resistance thermometer. The signal is read by the program through the temperature translator, the Amp, Phase, Temp card, the Relay card, and voltage monitor to the calculator. The measured temperature is compared to values in the calculator program which control power to the environmental chamber heaters. Power is supplied at levels sufficient to maintain temperature increase of about $1^{\circ}\text{C}/\text{min}$.

BLOCK DIAGRAM OF SYSTEM




Figure 2

The Tension Amplifier, Amp, Phase, Temp. card, Relay card, Motor Power Supply and Amplifier, and Voltage Monitor are all contained in the Multiprogrammer.

The following data values are measured and recorded on a magnetic tape cassette for eventual data reduction:

Strain phase angle, degrees
Strain amplitude, cm.
Sample length, cm.
Average static tension, dynes
Temperature °C
Stress phase angle, degrees
Stress amplitude, dynes

MEASUREMENT FEATURES

1. Temperature

Temperature is detected by a 100-ohm platinum resistance thermometer (RDF #80RB) and is measured and transmitted through the multiprogrammer to the calculator by a temperature transmitter card (RDF Model 2050). This card provides a linearized output signal of 1 mv per degree. Accuracy is $\pm 1^\circ\text{C}$ over a full scale range (-180°C to $+650^\circ\text{C}$).

The temperature chamber is designed with sufficient thermal mass to make control unnecessary from -150°C to approximately -60°C . Natural heat absorption after removal of liquid N_2 cooling causes the system to increase in temperature at a rate very close to $1^\circ/\text{min}$. which is about ideal for data acquisition. Above -30°C some power must be supplied to the heaters to maintain the rate of temperature rise. This is handled by a heater control system consisting of two solid state relays which proportion power to the heaters in amounts depending on needs. Power requirements are divided into four approximate levels:

5%, 30%, 50% and 70%. Each level covers a specific range of temperature and the program automatically switches to the next higher power when each temperature level is reached. The system maintains a rate of temperature increase of approximately $1^{\circ}\text{C}/\text{min}$.

2. Modulus and Phase Measurement

A lock-in amplifier (Princeton Applied Research Model 5204) is used to measure in-phase and quadrature amplitudes of both the load and displacement signals and the phase relation of each signal with respect to a reference from the Rheovibron. The calculator program sends a command to the multiprogrammer to switch on the load signal to the lock-in amplifier which, after a program delay for settling, reads the in-phase and quadrature components with respect to the reference, and stores them in the calculator memory.




Figure 3

After these values are stored, the multiprogrammer receives a command to switch the dynamic displacement signal into the lock-in amplifier where it repeats the same measurements on these signals and stores the results in the calculator:




Figure 4

The calculator program then computes the voltage magnitudes of both the load and displacement signals, P and X, and the phase angle, δ , between the load and displacement.

$$P = \sqrt{P_I^2 + P_Q^2}$$

$$X = \sqrt{X_I^2 + X_Q^2}$$

$$\delta = \phi_1 - \phi_2 = \tan^{-1} \frac{P_Q}{P_I} - \tan^{-1} \frac{X_Q}{X_I}$$

The load and displacement signals are scaled by the calibration factors stored in the program and the complex modulus, E, is computed and displayed on the calculator. At the same time the tangent of the loss angle, $\tan \delta$, is displayed with the temperature. This modulus is uncorrected for instrument compliance but is corrected for changes in sample length which will be discussed in the section on tensioning the sample.

3. Sample Tension Control

The type of sample most frequently used (thin films, fibers, monofilaments, etc.) on the Rheovibron will not support a compressive load without buckling. This condition would, of course, distort the negative portion of the sine wave making measurements impossible. To alleviate this situation the sample must be maintained under sufficient tension to eliminate buckling but insufficient to introduce non-linearity or excessive creep into the experiment. The amount of tension is critical to the measurement.

The tension control system uses a stepping motor (Superior Electric Slo-Syn Model 5850-1008) to add or subtract small increments of tension to the sample. The motor receives its commands from the calculator. Each step of the motor represents a change in length of the sample of 1.25×10^{-4} cm. The program computes the modulus of the sample and uses this value to establish maximum and minimum limits of tension. It then calculates the magnitude of a single step to determine its load value and computes the number of steps required to put the sample tension within limits. The program then instructs the motor, through the multiprogrammer, to initiate the calculated number of steps, either + or -. The tension is measured and corrected twice during the taking of each datum point. As the sample softens, the program changes tension limits. An up and down count is stored in a register which, when multiplied by the value of each step, is added to, or subtracted from the initial sample length.

TYPICAL TEST PROCEDURE

Following is a step by step description of a typical test.

1. Turn on power to all systems and allow at least a 20 minute warmup.
2. Prepare the sample using care that its size does not exceed the instrument force capability as explained on pages 41-44 of Rheovibron Manual #68. Thickness and width, in cm, are carefully measured and recorded.
3. Install sample in the grips using extreme care that it is straight and will not develop side thrusts or distortions when tension is applied. This is very important and has great bearing on quality of the data.
4. Load the control program into the calculator and strike RUN.
5. Type in the parameters requested by the program. After each entry, strike CONTINUE.
6. At the end of the parameters routine the calculator will display "Ready for Cooldown?". If yes, strike CONTINUE.
7. Add liquid nitrogen to the chamber and to the Dewar cooling the N₂ gas. Start gas flow at high rate, about 16-18 SCFH on the flowmeter. Keep a good supply of liquid N₂ in both the chamber container and the Dewar during cooldown.
8. The program will follow the course of cooldown, automatically adjusting tension as the sample shrinks. When starting temperature is reached, the calculator will display ...Minimum Temperature..., and the program will halt.
9. Reduce the nitrogen gas flow to about 6-8 SCFH and allow the liquid N₂ to boil off the temperature chamber tank. The main data taking program is initiated by striking CONTINUE.
10. From this point data taking and tension control are automatic. The temperature will rise approximately 1°C/min. from the heat obtained from surrounding. At about -60°C, remove the copper coil from the liquid N₂ Dewar. Power to the heaters for continuing the temperature rise of about 1°C/min. is handled by the program.
11. When any preset limit is reached, Max. Temperature, Max. Sample Length, or Max. Static Tension, the system will automatically halt.

12. During the course of the experiment, each measurement will be displayed in terms of E, $\tan \delta$ and Temperature for convenience of the operator.
13. After a run, data are handled by loading the Data Reduction Program into the calculator and allowing it to search the recorded tape for data.

CONCLUSIONS AND RECOMMENDATIONS

A working prototype system has been completed which achieved the following improvements:

<u>Parameter</u>	<u>Manual Operation</u>	<u>Automatic Operation</u>
Labor/hours/run	3 - 6	1/2 - 1
Resolution of $\tan \delta$ measurement	~.003	.001
Tan δ range	.001 - 1.7	.001 - >3
Time for data reduction	3 - 5 hours	a few minutes (2 - 10)

An attempt was made to keep "hardware" at a minimum and take advantage of the flexibility of the calculator-computer which is part of the system. The program as presently written works quite well but will be improved and refined with experience.

Possible areas for future work and development might be:

- a) Introduce servo-control (closed loop) on stress and/or strain amplitudes. This could extend the system to creep, stress relaxation and thermal expansion measurements.
- b) Write programs for digital printer/plotter.
- c) Develop hardware and write programs for automatic analog plotting of data.

COOLDOWN ROUTINE

This program is designed to maintain sufficient tension on the sample to keep it straight as it is cooled. The tension control senses buildup of force as the sample shrinks and the stepping motor steps to relieve tension to keep within the limits. The program constantly checks for excesses in any of the preset test limits and for the preset minimum temperature. When the minimum or starting temperature is reached, the display shows ...Minimum Temperature... and the program halts. To start the test run, strike CONTINUE.

COOLDOWN PROGRAM




Figure 5

MAIN CONTROL ROUTINE

The main control program is designed to measure and record data as the temperature is increased. It is a repetitive process which starts with a check on the size of the signal from the load transducer for the purpose of autoranging the lock-in amplifier. After autoranging the system reads the in-phase and quadrature load signals and stores them in the calculator memory. The program then branches into a subroutine to measure static tension on the sample and make any corrections, + or -, to keep it within proper limits. It then proceeds to measure and store the in-phase and quadrature displacement signals from the Rheovibron. Next, the program makes a correction in sample length for any + or - steps required to correct the tension. That is followed by another tension measurement and correction, if necessary. At this point the temperature is measured and stores. The program then makes use of the scale and calibration factors prestored in the memory to calculate the modulus, E , and $\tan \delta$. The values are displayed with temperatures so that the operator can follow the course of the test. The system then checks to see if limits have been exceeded and, if not, records the data on tape. The process then repeats.




Figure 6

DATA REDUCTION ROUTINE

The data reduction program utilizes a standard curve fitting subroutine that smooths raw data by means of a 4th order polynomial regression. The first step in the program is to ask the operator to select a file number that specifies the header for the experimental data to be analyzed. The program will check to see if this file is a valid header file. This is accomplished by identifying the file and testing to see if the actual byte size is 160 -- the size of a header file. If not, a message is displayed and the program halts. If it is a valid header file, the processing continues.

The next step is the calculation of the sample volume from the input specifications. If the volume is zero, the program halts.

Next the program sets the initial temperature reference and the temperature limit which is the temperature reference plus ten degrees. This allows for smoothing over a ten degree range.




Figure 7A

The program next reads the data sets, each consisting of a 7 element array containing the measured parameters needed to plot a tan delta curve. Young's Modulus and the angle, δ , are calculated from these data points and, along with the temperature, they are stored in arrays. When the temperature reading has exceeded the temperature limit for this curve the arrays are read into a subroutine "FIT2" which will calculate the compliance coefficient for the polynomial regression. These coefficients are then used in the equations to calculate "E" and tan delta. Before the next set of data prints are read from the tape the last five data sets in the arrays are stored into the first five locations of the next data sets. This will allow for a smooth progression. This process repeats itself for every ten degrees until the data runs out or until the maximum temperature has been reached.





Figure 7B


PARTS LISTS

DIAGRAMS

4P 9825A

H.P. 6340 P





Overall Diagram of Interconnections.

MAIN FRAME

Ref.Desig.	Description	Q	Manufacturer
			<u>Hewlett Packard:</u>
CALC	Calculator	1	H.P.#9825A
ROM	9872 Plotter-General I/O Extended I/O	1	H.P.#98216A
MUL	Multiprogrammer	1	H.P.#6940B
	Voltage Regulator	1	H.P.#69351B
RELAY	Relay Output/Readback Card	1	H.P.#69433A
	Voltage Monitor Card	1	H.P.#69421A
PAR	Lock-In Amplifier	1	Princeton Applied Research #5204



CABINET

Ref.Desig.	Description	Q	Manufacturer
	Frame, Base & Casters	1	EQUIPTO 40-3532-03
	Side Panel	2	EQUIPTO 30-3538-11
	Top Panel	1	EQUIPTO 20-2112-46
	Rear Panel, Louvered	1	EQUIPTO 20-2812-14
	Equipment Shelves	2	EQUIPTO 18-1925-01

MOTOR DRIVE ASSEMBLY

Ref.Desig.	Description	Q	Manufacturer
R1	Stepping Motor	1	Superior Elec. SS50-1008
	Resistor, 5Ω, 25W	1	Superior Elec. BP262-G10
	Connector	1	Amphenol 3102A-14S-06P
	<u>Mechanical Parts:</u>		
	Housing		IDR
	Front Plate		IDR
	Shaft, Precision	1	PIC A3-30
	Bearing	2	PIC B10-9
	Collars	2	PIC C1-3
	Pulley	1	PIC FD6-64
	Pulley	1	PIC FD6-18
	Hub Clamp	4	PIC L4-6
	Drive Belt	1	PIC FA-95
	Motor Mount Plate	1	IDR
	Spacers - Housing	2	IDR
	Spacer - Motor Shaft	1	IDR
	Bellows Coupling	1	PIC T1-10
	Shaft, Precision	1	PIC A3-50
	1" Louvere	1	Midget Louvere Co.


MOTOR DRIVE ASSEMBLY



-Pic L4-6


SLO-SYN DRIVE ASSEMBLY

MOTOR DRIVE ASSEMBLY



MOTOR DRIVE ASSEMBLY


Motor



CABLES

Ref.Desig.	Description	Q	Manufacturer
GPI/O	16-bit Duplex Interface	1	Hewlett Packard 98032A
RS232-C	Serial I/O Interface	1	with opt. 040 Hewlett Packard 98036A
Dyn. Load	2 Conductor-Shielded	1	IDR #101
Dyn. Load- Displ.	3 Conductor-Shielded	1	IDR #102
Dyn. Dyn. Load, Load/ Disp., Ref.	BNC-BNC 6 ft.	3	IDR #103
I, Q	BNC-open 6 ft.	2	IDR #104
Range	5 Conductor Cable - Shielded	1	IDR #105
RTD	Sheathed Temp. Probe	1	RDF PN21701/110
Motor	6 Conductor Shielded	1	IDR #106
Power Strip	6 Outlet Wired Plug Mold	1	Equipto 03-0606-32


Interconnecting Harness




AMP.-PHASE-TEMP. Card

Ref.Desig.	Description	Q	Manufacturer
A1	Amplifier	1	Analog Devices 183J
A2, A3	Amplifier	2	Analog Devices 741H
R1	Trim Pot, 50K Ω	1	Allen Bradley RT5
R2	Resistor, Metal Film, 1.33K Ω	1	Corning C4
R3, R5	Resistor, Metal Film, 1K Ω	2	Corning C4
R4	Resistor, Metal Film, 200K Ω	1	Corning C4
R6	Resistor, Metal Film, 4.42K Ω	1	Corning C4
R7	Resistor, Metal Film	1	Analog Devices - Matched
R8, R10	Resistor, Metal Film, 2K Ω	2	Corning C4
R9, R11	Resistor, Metal Film, 1.2K Ω	2	Corning C4
R12, R14	Resistor, Metal Film, 10K Ω	2	Corning C4
R13, R15, R16, R17	Trim Pot, 10K Ω	4	Allen Bradley RT5
R18, R19	Resistor, Metal Film, 60K Ω	2	Corning C4
R20, R21	Resistor, Metal Film, 10K Ω	2	Corning C4
R22, R23	Resistor, Metal Film, 51K Ω	2	Corning C4
C1, C2	Capacitor, 5.0 MFD, 50VDC, $\pm 5\%$	2	S&EI 22U
C3	Capacitor, 0.68 MFD, 50VDC, $\pm 5\%$	1	S&EI 22U
<u>Mechanical Parts:</u>			
	Connector 2x15	1	Hewlett Packard 251-15-30-261
	IC Socket 30 pin	2	Augat 8058-39G4
	Connector Jack	7	Cambion 450-3704
	P.C. Solder Terminal	50	Cambion 1461-2
	P.C. Board	1	IDR AV-100

AMP.-PHASE-TEMP. Card




AMP.-PHASE-TEMP. Card




TENSION CARD

Ref.Design.	Description	Q	Manufacturer
A1	Instrumentation Amplifier	1	Analog Devices 610L
R2	Trim Pot, 100K Ω	1	Allen Bradley RT5
R1	Resistor 402 Ω 1%	1	Corning RN65D
	<u>Mechanical Parts:</u>		
	Connector 2x15 30 pin	1	Hewlett Packard 251-15-30-261
	P.C. Solder Terminals	13	Cambion 1461-2
	P.C. Board	1	IDR AV-101

TENSION Card




TENSION Card




MOTOR CARD

Ref.Desig.	Description	Q	Manufacturer
Motor	SLO-SYN Translator	1	Superior Elec. STM 1800
R1, R3	Resistor, Metal Film, 100K Ω	2	Corning C4
R2	Resistor, Metal Film, 47 Ω	1	
R4	Resistor, Metal Film, 20 Ω 10W	1	
C1	Capacitor 0.05 MFD \pm 5% 50VDC	1	S&EI 22U
C2	Capacitor 0.1 MFD \pm 5% 50VDC	1	S&EI 22U
	<u>Mechanical Parts:</u>		
	Connector 2x15 30 pin	1	Hewlett Packard 251-15-30-261
	Connector Jack	16	Cambion 450-3704

MOTOR Card






RANGE CARD

Ref.Desig.	Description	Q	Manufacturer
O. Bread-board	P.C. Board	1	Hewlett Packard 69380A
	<u>In & P.A.R. Lock-In Amplifier</u>		
	K1,K2,K3,K4 Relay, SPST, 5VDC	4	Sigma 191TE1AZ-5S
	R1,R2,R3,R4 Resistors, Metal Film, 10K Ω	4	Corning C5
	Q1,Q2,Q3,Q4 Transistors	4	Motorola 2N4124
	<u>Mechanical Parts:</u>		
	P.S. Solder Terminals	25	Cambion 1461-2
	Epoxy-Glass Board	1	
	Mounting Bracket	1	
	Spacers	2	
	Screws	2	


RANGE Card

THIS PAGE IS BEST QUALITY PRACTICABLE
FROM COPY FURNISHED TO DDC



RANGE Card


THIS PAGE IS BEST QUALITY PRACTICABLE
FROM COPY FURNISHED TO DDG




RHEOVIBRON RELAY BOARD

Ref.Desig.	Description	Q	Manufacturer
K1, K2	4.5" Long Vector Dip Plugboard	1	Vector Electronics 3662-5
	Relay DPST	2	Magnecraft W171- DIP25

RELAY Board




RELAY Board



SV

RELAY Board




RHEOVIBRON INTERFACE

Ref.Desig.	Description	Q	Manufacturer
	Resistor, Metal Film, 100K Ω	1	Corning C4
	Resistor, Metal Film, 1.5K Ω	1	Corning C4
	Resistor, Metal Film, 10 Ω	1	Corning C4
	Trim Pot 100K		Allen Bradley RT5
	Capacitor, 1MFD, 100V 5%	1	S&EI 22U
	<u>Mechanical Parts:</u>		
	Box 4" x 4" x 2"	1	Bud AU1083HG
	BNC Bulkhead Conn.	3	Amphenol 31-236
	Connector	1	Amphenol 3102A-10SL-03P
	Connector	1	Amphenol 3102A-12S-03P

RHEOVIBRON Interface

THIS PAGE IS BEST QUALITY PRACTICABLE
FROM COPY FURNISHED TO DDG



HEATING CONTROL

Ref.Desig.	Description	Q	Manufacturer
0. Bread-board	Breadboard P.C. Board	1	Hewlett Packard 69380A
K1, K2	Relay SPST	2	Magnecraft W107DIP5 or Sigma 191TE1AZ-5S
Q1, Q2	Transistors 2N3903	2	Motorola M2N3903
R1, R2	Resistor, Metal Film, 10K Ω	2	Corning C4


**THIS PAGE IS BEST QUALITY PRACTICABLE
FROM COPY FURNISHED TO DDG**



HEATING CONTROL

THIS PAGE IS BEST QUALITY PRACTICABLE
FROM COPY FURNISHED TO DDG

98380 A




Sigma or Denver
relays

POWER SUPPLY PANEL

Ref.Design.	Description	Q	Manufacturer
K1, K2	Power Supply	1	Lamda LOS-Y-15
	Power Relay	2	Sigma 221A-1-5D
	Temperature Transmitter	1	RDF 2050-80RB
	<u>Mechanical Parts:</u>		
	Rear Panel	1	IDR #

POWER SUPPLY PANEL


THIS PAGE IS BEST QUALITY PRACTICABLE
FROM COPY FURNISHED TO DDG



Mounting for Temperature Chamber



Ref.Desig.	Description	Q	Manufacturer
	Milling Attachment	1	Clausing M6-500
	Bottom Plate	1	IDR AV-A125
	Adapters	2	IDR AV-A126

Mounting for Temperature Chamber



ONE REQUIRED
STEEL


THIS PAGE IS BEST QUALITY PRACTICABLE
FROM COPY FURNISHED TO DDC



ONE REQUIRED

Mounting for Temperature Chamber

(3) $\frac{5}{16}$ DIA



THIS PAGE IS BEST QUALITY PRACTICABLE
FROM COPY FURNISHED TO DDG

BASE

ONE RECEIVED


ALUMINUM

6061-T6

2024-T6

REFERENCES

1. Rheovibron Instruction Manual No. 62 (Nov. 1973)
2. Rheovibron Instruction Manual No. 68 (April 1973)
3. Rheovibron Maintenance Manual No. 68(70)-2 (June 1975)
4. Hewlett Packard 9825A Calculator Operating and Programming No. 09825-90000
5. Hewlett Packard 6940B Multiprogrammer Users Guide No. 59500-90005
6. Hewlett Packard Operating and Service Manual for Voltage Regulator Card, No. 69351-90002
7. Hewlett Packard Operating and Service Manual for Relay Output/Readback Card, No. 69433-90001
8. Hewlett Packard Operating and Service Manual for Voltage Monitor Card, No. 69421-90001
9. Hewlett Packard Operating and Service Manual for Breadboard Output Card, No. 69380-90001
10. Hewlett Packard Extended I/O Programming, No. 09825-90025
11. Hewlett Packard Installation and Service Manual for 16-Bit Interface, No. 98032-900006
12. Hewlett Packard Operating and Service Manual for Multiprogrammer, No. 06940-9000
13. Operating and Service Manual for Princeton Applied Research Lab. Lock-In Amplifier Model 5204, No. M5204; 8/77-GG
14. Hewlett Packard Installation and Service Manual for Serial I/O Interface, No. 98036-9-9000



PROCEDURES MANUAL

AUTOMATED RHEOVIBRON

CALIBRATION PROCEDURE

Using instructions contained in Rheovibron Manual #68, pages 14-18, balance, zero and calibrate both the strain and the stress gages.

Set up the instrument controls as follows:


- Bal-Calib Switch on Meas
- Main Switch on AF
- Tan δ Range Switch on 30
- AF Switch on 30
- DF Dial on 1000
- Amp. Adjust Dial fully counterclockwise
- Phase and Meas Switch on Meas
- Frequency Selector on any frequency above 1 Hz

Turn power on to the calculator and the multiprogrammer and load the program into the calculator memory. This is accomplished by inserting the mag tape cartridge, keying in fdf ϕ , striking the EXECUTE key, and then keying in and executing the command to load fdf ϕ . A pilot lamp next to the tape cartridge will light as the program is being loaded into the calculator. Run the "ZEROES:" program by keying into calculator the command RUN 210, striking the EXECUTE key. When the calculator requests a Z, key in number 6, Continue. A number will appear on the display which will represent in-phase component of the output of the strain gage (X_I). This should be zero (± 0.005 volts). If it is not, refer to the instructions for removing offsets from the system.

When the in-phase channel is satisfactorily zeroed, repeat the program by keying RUN 210 into the calculator and striking EXECUTE. When a Z is requested, key in 12, Continue. You will now read the quadrature component of the strain (X_Q) and it, too, should be zero (± 0.005 volts).

With no sample in the grips, gradually turn the Amp. Adjust knob clockwise until the tan δ meter registers exactly full scale. Again run the "ZEROES:" program by keying Run 210, EXECUTE. Key in a 6 for Z.

The number which appears on the display represents the magnitude of the in-phase strain component, X_I . Record it. Next rerun the program by keying RUN 210, EXECUTE, and providing a 12 for the Z. The number which now appears on the display represents the quadrature component of the strain transducer, X_Q . From these two values, X_I and X_Q , we can calculate the magnitude of the output of the strain transducer under these conditions:




$$|X| = \sqrt{X_I^2 + X_Q^2}, \text{ volts}$$

Figure 8

The actual dynamic displacement amplitude measured by the transducer is exactly 5.0×10^{-3} cm. The calculation is explained on page 6 of the Rheovibron Manual #68. $|X| \div 5.0 \times 10^{-3} =$ scale factor for strain transducer, V./cm. This is the strain transducer calibration which is stored in the program.

Use the same procedure for the stress transducer except that a sample must be placed in the grips and sufficient tension applied to eliminate buckling of the sample when dynamic displacement is applied. The sample material is not important. The instrument controls are set the same as the previous measurement except that the main switch is now on DF instead of AF.

Gradually turn the Amp. Adjust control clockwise until the $\tan \delta$ meter reads exactly full scale. The in-phase component can now be read by keying in and executing RUN which will initiate the "ZEROES:" program. Provide a 5 for the Z which will display the in-phase component of load, P_I . The quadrature component is determined by again running the program and providing an 11 for Z. The calculation is the same:



$$|P| = \sqrt{P_I^2 + P_Q^2}, \text{ volts}$$

Figure 9

The load signal is now exactly 10^4 dynes and the determination is again explained in Rheovibron Manual #68, page 6.

$$\frac{|P|}{10^4} = \text{scale factor for the force transducer, V./dyne.}$$

This is the stress transducer calibration factor which is stored in the program.

PROCEDURE FOR CHECKING AND ADJUSTING SYSTEM ZEROES

1. Turn on the system.
2. Turn Amplitude Adjust fully counterclockwise.
3. Load the control program into the calculator.
4. Run the "Zeroes:" subroutine, RUN, EXECUTE.
5. When the "Z" is requested by the subroutine, enter 5 and strike "continue." This command will cause the system to switch on the stress signal and read the in-phase component P_I . The signal should be 0.000 ± 0.005 volts. Record the value.
6. Strike RESET, CONTINUE, which will cause the system to request another Z. Type in 6 and strike CONTINUE. The output of the in-phase strain component, X_I , will be displayed. It, too, should be 0.000 ± 0.005 volts. Record the value.
7. Repeat the same procedure using


<u>"Z"</u>	<u>channel displayed</u>
11	quadrature stress, P_Q
12	quadrature strain, X_Q

These values should be 0.000 ± 0.005 volts.

8. If any channel shows values out of the specified range, check the zero offset for that channel on the front of the lock-in amplifier. The + or - switch for zero offset must be depressed depending on the sign of the offset to be removed. Offsets are minimized from the amplifier by rotating the zero-offset potentiometer for that channel, searching for a minimum.
9. If the offset cannot be removed with lock-in amplifier controls, it will be necessary to trim the "zero" potentiometers on the Temperature, Phase, Amplitude board in the multiprogrammer. Refer to diagram on page 25 - R20 for in-phase adjustments and R21 for quadrature. You will find these trim pots on the board marked "Z". The trim pot on the upper edge of the board is for the in-phase adjustment, and the one on the lower edge of the board is for the quadrature adjustment.

CARD LOCATIONS
MULTIPROGRAMMER

<u>Card</u>	<u>Location</u>
Std. Input	100
Remote/Local	200
Logic & Timing	300
Relay	402B
Voltage Monitor	403C
Amplitude/Phase/Temp.	405E
Tension	407G
Heater	409I
Motor	411K
Range	414N
Unit Select	500
Voltage Reg.	600



CONTROL PROGRAM

Comments

accepting inputs of parameters from the operator

```
0: dsb "RHEOVIBR
ON DATA ACQUISSI
TION"
1: esb "PARAMETE
RS"
2: esb "COOLDOWN"
3: esb "GO"
4: "END":dsb
  "Not a null
  file - restart"
  iend
5: "PARAMETERS":
  dim R[20]
6: ent "File #
  for this test
  ?",r24
7: trk 0;fdf
  r24;idf r32;
  r33;if r33#0;
  sto "END"
8: mrk 1;160
9: ent "Sample
  #",r31;r31+R[1]
10: ent "Day",
  r33;r33+R[2]
11: ent "Month",
  r32;r32+R[3]
12: ent "Year",
  r34;r34+R[4]
13: ent "Sample
  length",r35;
  r35+R[5]
14: ent "Sample
  width",r36;r36+
  R[6]
15: ent "Sample
  thickness",r37;
  r37+R[7]
16: ent "Max.
  temp.",r16;r16+
  R[8]
17: ent "Min.
  temp.",r17;r17+
  R[9]
18: 3.797e3+r10;
  r10+R[10]
19: 6.53e-4+r11;
  r11+R[11]
20: ent "Frequen
  cy",r38;r38+R[1
  2]
21: ent "C1",
  r41;r41+R[13]
22: ent "C2",
  r42;r42+R[14]
23: ent "C3",
  r43;r43+R[15]
24: ent "Max. loa
  d",r18;r18+R[16]
]
```

Routine for selecting tape file no. for recording data and accepting test parameters.


```

25: ent "Max.
    sample length",
    r15;r15+R[17]
26: ent "Min.
    stress amplitud
    e",r14;r14+R[18]
    ]
27: fdr r24;rcf
    r24,R[*]
28: fmt c5,f9.0
29: wrt 16,"smp#
    ",r31
30: spc
31: fmt c3,2x,
    f2.0,c,f2.0,c,
    f4.0
32: wrt 16,"dts"
    ,r32,"/",r33,"/
    ",r34
33: spc
34: fmt c4,e9.3,
    c3
35: wrt 16,"leh=
    ",r35,"cm"
36: wrt 16,"udt=
    ",r36,"cm"
37: wrt 16,"thc=
    ",r37,"cm"
38: fmt c3,e9.2,
    c4
39: wrt 16,"+
    T=",r16,"C"
40: wrt 16,"-
    T=",r17,"C"
41: fmt c4,e9.3,
    c3
42: wrt 16,"coP=
    ",r10,"d/V"
43: wrt 16,"coX=
    ",r11,"c/V"
44: wrt 16,"fra=
    ",r38,"Hz"
45: wrt 16,"C1="
    ,r41
46: wrt 16,"C2="
    ,r42
47: wrt 16,"C3="
    ,r43
48: wrt 16,"mxP=
    ",r18,"a"
49: wrt 16,"mxL=
    ",r15,"cm"
50: wrt 16,"MIP=
    ",r14,"a"
51: spc 3
52: dsp "READY
    FOR COOLDOWN???
    ???"
53: stp
54: ret

```

Comments

These lines establish printer tape formats for printing the parameters entered in earlier steps.

```

55: "COOLDOWN":
56: le10+r8
57: esb "PC4"
58: esb "RT"
59: if r17>r5:
    jmp 3
60: esb "STP"
61: jmp -5
62: dsp "***Mini
    mum temperature
    ***"
63: stp
64: ret
65: "GO":
66: flt 3
67: "ALP":-200+r
    17:1+W
68: esb "RNG"
69: moct:wtb 9,
    170160,20001
70: dsp "I and
    Q stress readin
    g"
71: wait 15000
72: esb "PC4"
73: wait 15000
74: moct:wtb 9,
    170160,20005
75: .005+D
76: esb "R"
77: A+r27
78: moct:wtb 9,
    170160,20011
79: .005+D
80: esb "R"
81: A+r26
82: 2+W
83: moct:wtb 9,
    170040,160001
84: moct:wtb 9,
    170160,20006
85: dsp "I and
    Q strain readin
    g"
86: wait 30000
87: .005+D
88: esb "R"
89: A+r22
90: moct:wtb 9,
    170160,20012
91: .005+D
92: esb "R"
93: A+r21
94: r35+1.25e-4*
    r25+r3
95: (F+E)/2+r4
96: esb "PC4"
97: esb "RT"
98: esb "TC0"
99: esb "IODE"
100: esb "PRT"
101: esb "STP"
102: esb "RCF"
103: eto "ALP"
104: ret

```

Comments

Routine for controlling sample tension and measuring temperature during cooldown in preparation for a test run.

The main routine for collecting data

Reading the load transducer

Reading the displacement transducer

94: Correcting sample length to be used in the modulus calculation

96: Correcting sample tension

97: Reading temperature

98: Adjusting power to heaters

99: Calculating E and tan δ

100: Printing values of Temperature, E and tan δ

101: Checking for limits

102: Recording data

```

105: "RT":
106: dsb "Temper
ature reading"
107: 0+N*B
108: moct!wtb 9,
170160,20020+M
109: -.155+D;
asb "R"
110: N+1+N;A+
5+B;B/N+r5
111: if N<100;
jmp -2
112: ret
113: "R":
114: 0+A
115: moct!wtb 9,
170260,30000,
170240
116: wti 0,11;
wti 4,30000;
rdi 4+A;mdac
117: otdA+A;if
A>2047;A-4096+A
118: D*A+A
119: ret
120: "RDF":
121: r24+1+r24
122: trk 0;fdf
r24;mrk 1,56
123: rcf r24,r1,
r7
124: ret
125: "PM":
126: dsb "Tensio
n reading"
127: 0+N;100+E;
.005+D;0+F
128: moct!wtb 9,
170160,20040+M
129: asb "R"
130: if A<E;A+E
131: if A>F;A+F
132: N+1+N;if
N<160;jmp -3
133: ret
134: "RNG":
135: asb "PM"
136: F-E+G;G/
2+G;.707G+G
137: if G>.9;
25+r20;160010+C
;jmp 4
138: if G>.23;
10+r20;160004+C
;jmp 3
139: if G>.09;
2.5+r20;160002+
C;jmp 2
140: if G<.09;
1+r20;160001+C
141: moct!wtb 9,
170040,C
142: ret

```

Comments

Routine for measuring temperature

Routine for reading the multiprogrammer

Routine for recording data

Measuring minimum and maximum of
load sine wave for establishing
tension limits

Routine for autoranging the
lock-in amplifier

Comments

```
143: "PRI":  
144: flt 3  
145: prt "T=";r5  
146: prt "E";r8  
147: prt "tan";  
    r9  
148: spc 3  
149: ret  
150: "STP":if  
    r15<r3;prt "max  
    .L limit";jmp 5  
151: if r16<r5;  
    prt "max.temp.1  
    imit";jmp 4  
152: if r17>r5;  
    prt "min.temp.1  
    imit";jmp 3  
153: if r4>r18;  
    prt "max.load  
    limit";jmp 2  
154: jmp 5  
155: moctiwtb 9,  
    170000,27777  
156: fmt 7x;c13  
157: wrt 0,"***  
    LIMITS***"  
158: stp  
159: ret  
160: "+S":  
161: moctiwtb 9,  
    170160,20201,  
    21001  
162: wtb 9,17004  
    0,20000+W  
163: r25+1+r25  
164: ret  
165: "-S":  
166: moctiwtb 9,  
    170160,20200+W,  
    20400+W  
167: wtb 9,17004  
    0,20000+W  
168: r25-1+r25  
169: ret
```

Routine for printing T, E, and tan δ

Routine for checking limits. Program
halts if limits are exceeded.

Routine for stepping the motor to
add to sample tension.

Routine for subtracting from sample
tension.

```

170: "IQDE":
171: r26/r27+A+r
    6:atn(r6)+r6
172: r21/r22+B+r
    1:atn(r1)+r1
173: atn(A)-atn(
    B)+r19:abs(r19)
    +r19
174: tan(r19)+r3
175: flt 3
176: r(r27+2+
    r26+2)+r28
177: r28+r20+
    r10+r7
178: r(r22+2+
    r21+2)+r23
179: r23+r11+r2
180: r7/r2*r3+2+
    A
181: r35*r36+
    r37+B
182: A/B+r8
183: 1.24e-4/r3+
    r8*(r36+r37)+Y+
    2.09e-5*Y+Y
184: Y+r2+Y+.5Y+
    Y
185: flt 3
186: ret
187: "PC4":
188: 0+N
189: "DLP":N+1+N
190: if N>3:sto
    "CLP"
191: r8*2.4432e-
    12+Y
192: r8*1e-10+r1
    3
193: r13+10Y+r12
194: esb "PM"
195: if E<r13:
    jmp 3
196: if E>r12:
    jmp 7
197: sto "BLP"
198: (r13-E)/
    Y+B:abs(B)+1+B:
    fxd 0
199: B-1+B
200: esb "+S"
201: if B>0:jmp
    -2
202: sto "DLP"
203: if Y<1e-2:
    sto "CLP"
204: (E-r12)/
    Y+B:abs(B)+1+B:
    fxd 0
205: B-1+B
206: esb "-S"
207: if B>0:jmp
    -2
208: "BLP":flt
    3:sto "DLP"
209: "CLP":ret

```

Comments

Routine for calculating E
and tan δ

Routine for adjusting sample
tension

```

210: "ZEROS":
211: ent Z
212: moct:wtb 9,
    170160,20000+Z
213: .005+0:esb
    "R"
214: fxd 3:dsb R
215: ato 2
216: "CHECK1":
217: ent Y
218: moct:wtb 9,
    170040,16000+Y
219: atp
220: moct:wtb 9,
    170160,16000
221: ato 10
222: end
223: "AMPS":
224: ent Z
225: moct:wtb 9,
    170160,20000+Z,
    30000
226: moct:wtb 9,
    170040,160002
227: .005+0:esb
    "R"
228: fxd 3:dsb R
229: ato 227
230: "TC0":
231: if r5>200:
    moct:wtb 9,1700
    40,110024;jmp
    13
232: if r5>120:
    moct:wtb 9,1700
    40,110004;jmp
    12
233: if r5>60:
    jmp 3
234: if r5>20:
    jmp 5
235: if r5>-30:
    jmp 6
236: K-1+K
237: if K=1:moct
    :wtb 9,170040,
    110000;jmp 2
238: if K=0:moct
    :wtb 9,170040,
    110004;2+K;jmp
    6
239: L-1+L
240: if L=1:moct
    :wtb 9,170040,
    110004;2+L;jmp
    4
241: moct:wtb 9,
    170040,110004
242: wait 10000
243: moct:wtb 9,
    170040,110000
244: ret
*13633

```

Comments

Routine for checking system zero

This routine is not used.

Routine for checking zero and calibration

Routine for adjusting power to the heaters in temperature chamber

RHEOVIBRON PROGRAM r Variables

r1 degrees, strain phase angle
r2 cm, strain amplitude
r3 cm, sample length
r4 dynes, average tension static
r5 °C, temperature
r6 degrees, stress phase angle
r7 dynes, stress amplitude
r8 dynes/cm², E
r9 tan δ
r10 dynes/volt, load calibration
r11 cm/volt, displacement calibration
r12 dynes, PC Max. load-tension control
r13 dynes, PC Min. load-tension control
r14 dynes, Min. stress amplitude
r15 cm, Max. sample length
r16 °C, Max. temperature
r17 °C, Min. temperature
r18 dynes, Max. load
r19 degrees, δ
r20 ranging scale factor "RNG":
r21 volts, quadrature strain, X_Q
r22 volts, in-phase strain, X_T
r23 volts, strain amplitude, X
r24 recording file register "RCF":
r25 "+S" & "-S" counter

r26 volts, quadrature stress, P_Q
r27 volts, in-phase stress, P_I
r28 volts, stress amplitude, P_V
r29
r30
r31 sample #
r32 month
r33 day
r34 year
r35 cm, sample length (initial)
r36 cm, sample width
r37 cm, sample thickness
r38 Hz, frequency
r39
r40
r41 C1
r42 C2
r43 C3